


# ■ DOA SERIES ■

## Urine Drug Test



Test Strip

Test Card

Single-Drug Urine Test



Premium Cup

Square Cup

Test Card

Select Cup

Hex Cup

Multi-Drug Urine Test

## Saliva Drug Test



Single-Drug Saliva Test Midstream



D-Pen

Test Card

Test Cup

Lollipop Test

Multi-Drug Saliva Test

# ■ DOA SERIES ■

## DOA Reader



### DC908A

- Portable and easy to carry
- Barcode verification
- Capture results in 30 seconds
- Work with PC interface through the connection of USB port

### DC911A


- Portable and easy to carry
- Rapid, multi-assay, on-the-spot analyzer
- Data management
- Equipped with a display screen for clearer results
- Built-in printer

### DC909A

- Portable and easy to carry
- More objective than visual interpretation
- Capture results in 30 seconds
- Storage more than 10000 datas

# ■ DOA SERIES ■

## Hair Drug Test



- Capture results in 5 minutes
- Long detection window for 90 days
- Quantitative detection of drugs in hair

## Test Items Available

Specimen	Test Item	Cut-off level	Format	Qualification
Urine	Amphetamine (AMP)	1000/500/300 ng/mL	Strip/Card	CE marked
	Barbiturates (BAR)	300 ng/mL	Strip/Card	CE marked
	Benzodiazepines (BZO)	300/200/100 ng/mL	Strip/Card	CE marked
	Buprenorphine (BUP)	10 ng/mL	Strip/Card	CE marked
	Cocaine (COC)	300/150/100 ng/mL	Strip/Card	CE marked
	Marijuana (THC)	50/40/25 ng/mL	Strip/Card	CE marked
	Methamphetamine (MET)	1000/500/300 ng/mL	Strip/Card	CE marked
	Methadone (MTD)	300/200 ng/mL	Strip/Card	CE marked
	Methylenedioxymethamphetamine - ecstasy (MDMA)	500/300 ng/mL	Strip/Card	CE marked
	Morphine (MOP)	300/100 ng/mL	Strip/Card	CE marked

# ■ DOA SERIES ■

## Test Items Available

Specimen	Test Item	Cut-off level	Format	Qualification
Urine	Oxycodone (OXY)	100 ng/mL	Strip/Card	CE marked
	Phencyclidine (PCP)	25 ng/mL	Strip/Card	CE marked
	Propoxyphene (PPX)	300 ng/mL	Strip/Card	CE marked
	Tri-cyclic Antidepressants (TCA)	1000 ng/mL	Strip/Card	CE marked
	Synthetic Cannabis (K2)	50/25 ng/mL	Strip/Card	CE marked
	Opiate (OPI)	2000 ng/mL	Strip/Card	CE marked
	Tramadol HCl (TRA)	1000/200 ng/mL	Strip/Card	CE marked
	Methadone Metabolite (EDDP)	300 ng/mL	Strip/Card	CE marked
	Fentanyl (FTY)	20 ng/mL	Strip/Card	CE marked
	Cotinine (COT)	200 ng/mL	Strip/Card	CE marked
	Ethyl Glucuronide (ETG)	500 ng/mL	Strip/Card	CE marked
	Ketamine (KET)	1000/500 ng/mL	Strip/Card	CE marked
	6-MAN (Monoacetylmorphine)	10 ng/mL	Strip/Card	CE marked
Saliva	Amphetamine (AMP)	50/40 ng/mL	Card/Midstream	CE marked
	Barbiturates (BAR)	60/20 ng/mL	Card/Midstream	CE marked
	Benzodiazepines (BZO)	30/10 ng/mL	Card/Midstream	CE marked
	Buprenorphine (BUP)	5 ng/mL	Card/Midstream	CE marked
	Cocaine (COC)	50/20 ng/mL	Card/Midstream	CE marked
	Marijuana (THC)	50/40/12 ng/mL	Card/Midstream	CE marked
	Methadone (MTD)	30 ng/mL	Card/Midstream	CE marked
	Methamphetamine (MET)	50 ng/mL	Card/Midstream	CE marked
	Methylenedioxymethamphetamine (MDMA)	100/50/25 ng/mL	Card/Midstream	CE marked
	Opiate (OPI)	50/40/15 ng/mL	Card/Midstream	CE marked
	Morphine (MOP)	15 ng/mL	Card/Midstream	CE marked
	Oxycodone (OXY)	20 ng/mL	Card/Midstream	CE marked
	Phencyclidine (PCP)	10 ng/mL	Card/Midstream	CE marked
	Cotinine (COT)	20/50 ng/mL	Card/Midstream	CE marked
	Propoxyphene (PPX)	25 ng/mL	Card/Midstream	CE marked
	Alcohol	0.02%	Card/Midstream	CE marked
	ETG	300/500 ng/mL	Card/Midstream	CE marked
Hair	Morphine (MOP)	0.2 ng/mg	Card	CE marked
	Methamphetamine (MET)	0.2 ng/mg	Card	CE marked
	Ketamine (KET)	0.2 ng/mg	Card	CE marked
Powder/Liquid	3 in 1 (MOP/MET/KET)	0.2 ng/mg	Card	CE marked
	Fentanyl (FTY)	10 ng/mL	Strip	CE marked
	Xylazine (XYL)	500 ng/mL	Strip	CE marked
	FTY+ XYL combo test	10ng/mL + 500 ng/mL	Strip	CE marked

*Note: More test items are coming.*